



danuButton® Compact XS / Compact Set / Premium Set

The danuButton® is a replacement feeding tube for long-term enteral nutrition and use with completely healed stoma. The low-profile feeding tube is easy to conceal under clothing and especially for active, mobile patients. The extensive accessories of the danuButton® Sets enable you to an easy and comfortable replacement of the gastrostomy probe, additional components are no longer necessary.

The decision for the danuButton® is also the decision for more quality of life through an innovative, high-quality product.



COMFORT

- **Large selection of sizes and lengths** for the individual and perfect fit
- **Extensive range of accessories** in the Sets for easy and convenient replacement of the gastrostomy tube – additional components not necessary
- **High comfort** due to soft, flexible silicone
- **Low-profile design** for more mobility and quality of life
- **Easy cleaning** and maintenance
- **Tilted** balloon valve

SAFETY


- Made of **high-quality medical grade silicone** with high clarity
- All sets are **latex-free, PVC-free, DEHP-free, BPA-free**
- **Radiopaque stripe** along the entire shaft
- **Apple-shaped inner retention balloon with high balloon stability** for long usage, secure fixation and to avoid leakage
- **Flat, rounded tip** enables atraumatic insertion
- **Recessed tip** to avoid trauma
- Secure **anti-reflux** valve

The danumed® Extension Sets included in the danuButton® Compact and Premium Sets are made of kink-resistant, biocompatible polyurethane. They are equipped with ENFit® connectors according to DIN EN ISO 80369-3 and danumed's innovative rotating connector, the ENSwivel®.



The danuButton® Compact Set contains:

- danuButton® Gastrostomy Button
- danumed® Extension Set for Gastrostomy Button, 30 cm, Y-ENFit®, ENSwivel®, angled button connector
- danumed® Introducer
- danumed® Balloon Syringe 10 ml
- optiLube lubricant
- Gauze Pads (4x)



| | | Fr-Size | | | | | |
|-----|--------------|--------------|--------------|--------------|--------------|--------------|----|
| | | 12 | 14 | 16 | 18 | 20 | 24 |
| 0.8 | DA 1042 1208 | DA 1042 1408 | DA 1042 1608 | DA 1042 1808 | DA 1042 2008 | DA 1042 2408 | |
| 1.0 | DA 1042 1210 | DA 1042 1410 | DA 1042 1610 | DA 1042 1810 | DA 1042 2010 | DA 1042 2410 | |
| 1.2 | DA 1042 1212 | DA 1042 1412 | DA 1042 1612 | DA 1042 1812 | DA 1042 2012 | DA 1042 2412 | |
| 1.5 | DA 1042 1215 | DA 1042 1415 | DA 1042 1615 | DA 1042 1815 | DA 1042 2015 | DA 1042 2415 | |
| 1.7 | DA 1042 1217 | DA 1042 1417 | DA 1042 1617 | DA 1042 1817 | DA 1042 2017 | DA 1042 2417 | |
| 2.0 | DA 1042 1220 | DA 1042 1420 | DA 1042 1620 | DA 1042 1820 | DA 1042 2020 | DA 1042 2420 | |
| 2.3 | DA 1042 1223 | DA 1042 1423 | DA 1042 1623 | DA 1042 1823 | DA 1042 2023 | DA 1042 2423 | |
| 2.5 | DA 1042 1225 | DA 1042 1425 | DA 1042 1625 | DA 1042 1825 | DA 1042 2025 | DA 1042 2425 | |
| 2.7 | DA 1042 1227 | DA 1042 1427 | DA 1042 1627 | DA 1042 1827 | DA 1042 2027 | DA 1042 2427 | |
| 3.0 | DA 1042 1230 | DA 1042 1430 | DA 1042 1630 | DA 1042 1830 | DA 1042 2030 | DA 1042 2430 | |
| 3.5 | DA 1042 1235 | DA 1042 1435 | DA 1042 1635 | DA 1042 1835 | DA 1042 2035 | DA 1042 2435 | |
| 4.0 | DA 1042 1240 | DA 1042 1440 | DA 1042 1640 | DA 1042 1840 | DA 1042 2040 | DA 1042 2440 | |
| 4.5 | DA 1042 1245 | DA 1042 1445 | DA 1042 1645 | DA 1042 1845 | DA 1042 2045 | DA 1042 2445 | |
| 5.0 | DA 1042 1250 | DA 1042 1450 | DA 1042 1650 | DA 1042 1850 | DA 1042 2050 | DA 1042 2450 | |
| 5.5 | | DA 1042 1455 | | | | | |
| 6.0 | | DA 1042 1460 | | DA 1042 1860 | | | |
| 6.5 | | DA 1042 1465 | | | | | |

1 Set per box, 5 boxes per shipping box.

The danuButton® Compact XS contains:


- danuButton® Gastrostomy Button
- danumed® Introducer
- danumed® Balloon Syringe 10 ml

| | | Fr-Size | | | | | |
|-----|--------------|--------------|--------------|--------------|--------------|--------------|----|
| | | 12 | 14 | 16 | 18 | 20 | 24 |
| 0.8 | DA 1043 1208 | DA 1043 1408 | DA 1043 1608 | DA 1043 1808 | DA 1043 2008 | DA 1043 2408 | |
| 1.0 | DA 1043 1210 | DA 1043 1410 | DA 1043 1610 | DA 1043 1810 | DA 1043 2010 | DA 1043 2410 | |
| 1.2 | DA 1043 1212 | DA 1043 1412 | DA 1043 1612 | DA 1043 1812 | DA 1043 2012 | DA 1043 2412 | |
| 1.5 | DA 1043 1215 | DA 1043 1415 | DA 1043 1615 | DA 1043 1815 | DA 1043 2015 | DA 1043 2415 | |
| 1.7 | DA 1043 1217 | DA 1043 1417 | DA 1043 1617 | DA 1043 1817 | DA 1043 2017 | DA 1043 2417 | |
| 2.0 | DA 1043 1220 | DA 1043 1420 | DA 1043 1620 | DA 1043 1820 | DA 1043 2020 | DA 1043 2420 | |
| 2.3 | DA 1043 1223 | DA 1043 1423 | DA 1043 1623 | DA 1043 1823 | DA 1043 2023 | DA 1043 2423 | |
| 2.5 | DA 1043 1225 | DA 1043 1425 | DA 1043 1625 | DA 1043 1825 | DA 1043 2025 | DA 1043 2425 | |
| 2.7 | DA 1043 1227 | DA 1043 1427 | DA 1043 1627 | DA 1043 1827 | DA 1043 2027 | DA 1043 2427 | |
| 3.0 | DA 1043 1230 | DA 1043 1430 | DA 1043 1630 | DA 1043 1830 | DA 1043 2030 | DA 1043 2430 | |
| 3.5 | DA 1043 1235 | DA 1043 1435 | DA 1043 1635 | DA 1043 1835 | DA 1043 2035 | DA 1043 2435 | |
| 4.0 | DA 1043 1240 | DA 1043 1440 | DA 1043 1640 | DA 1043 1840 | DA 1043 2040 | DA 1043 2440 | |
| 4.5 | DA 1043 1245 | DA 1043 1445 | DA 1043 1645 | DA 1043 1845 | DA 1043 2045 | DA 1043 2445 | |
| 5.0 | DA 1043 1250 | DA 1043 1450 | DA 1043 1650 | DA 1043 1850 | DA 1043 2050 | DA 1043 2450 | |
| 5.5 | | DA 1043 1455 | | | | | |
| 6.0 | | DA 1043 1460 | | DA 1043 1860 | | | |
| 6.5 | | DA 1043 1465 | | | | | |

1 Set per box, 5 boxes per shipping box.

The danuButton® Premium Set contains:

- danuButton® Gastrostomy Button
- danumed® Extension Set for Gastrostomy Button, 30 cm, Bolus-ENFit®, ENSwivel®, straight button connector
- danumed® Extension Set for Gastrostomy Button, 30 cm, Y-ENFit®, ENSwivel®, angled button connector
- danumed® Introducer
- danumed® Balloon Syringe 10 ml
- danumed® Enteral Syringe 60 ml, ENFit®
- optiLube lubricant
- Gauze Pads (4x)



| | | Fr-Size | | | | | |
|-----|--------------|--------------|--------------|--------------|--------------|--------------|----|
| | | 12 | 14 | 16 | 18 | 20 | 24 |
| 0.8 | DA 1041 1208 | DA 1041 1408 | DA 1041 1608 | DA 1041 1808 | DA 1041 2008 | DA 1041 2408 | |
| 1.0 | DA 1041 1210 | DA 1041 1410 | DA 1041 1610 | DA 1041 1810 | DA 1041 2010 | DA 1041 2410 | |
| 1.2 | DA 1041 1212 | DA 1041 1412 | DA 1041 1612 | DA 1041 1812 | DA 1041 2012 | DA 1041 2412 | |
| 1.5 | DA 1041 1215 | DA 1041 1415 | DA 1041 1615 | DA 1041 1815 | DA 1041 2015 | DA 1041 2415 | |
| 1.7 | DA 1041 1217 | DA 1041 1417 | DA 1041 1617 | DA 1041 1817 | DA 1041 2017 | DA 1041 2417 | |
| 2.0 | DA 1041 1220 | DA 1041 1420 | DA 1041 1620 | DA 1041 1820 | DA 1041 2020 | DA 1041 2420 | |
| 2.3 | DA 1041 1223 | DA 1041 1423 | DA 1041 1623 | DA 1041 1823 | DA 1041 2023 | DA 1041 2423 | |
| 2.5 | DA 1041 1225 | DA 1041 1425 | DA 1041 1625 | DA 1041 1825 | DA 1041 2025 | DA 1041 2425 | |
| 2.7 | DA 1041 1227 | DA 1041 1427 | DA 1041 1627 | DA 1041 1827 | DA 1041 2027 | DA 1041 2427 | |
| 3.0 | DA 1041 1230 | DA 1041 1430 | DA 1041 1630 | DA 1041 1830 | DA 1041 2030 | DA 1041 2430 | |
| 3.5 | DA 1041 1235 | DA 1041 1435 | DA 1041 1635 | DA 1041 1835 | DA 1041 2035 | DA 1041 2435 | |
| 4.0 | DA 1041 1240 | DA 1041 1440 | DA 1041 1640 | DA 1041 1840 | DA 1041 2040 | DA 1041 2440 | |
| 4.5 | DA 1041 1245 | DA 1041 1445 | DA 1041 1645 | DA 1041 1845 | DA 1041 2045 | DA 1041 2445 | |
| 5.0 | DA 1041 1250 | DA 1041 1450 | DA 1041 1650 | DA 1041 1850 | DA 1041 2050 | DA 1041 2450 | |
| 5.5 | | DA 1041 1455 | | | | | |
| 6.0 | | DA 1041 1460 | | DA 1041 1860 | | | |
| 6.5 | | DA 1041 1465 | | | | | |

1 Set per box, 5 boxes per shipping box.



All danuButton® Sets are:

- Latex-free
- PVC-free
- DEHP-free
- BPA-free