



# danuButton® Gastrostomy Button



**Safety meets  
Comfort**

German Manufacturer ✓

Reliable Supply Chain ✓

88 References ✓

Accessories in the Set ✓

High Patient Comfort ✓

Long-lasting ✓

Secure Hold ✓

Apple-shaped ✓

• Latex-free • PVC-free • DEHP-free • BPA-free

# danuButton® Gastrostomy Button


The danuButton® is a low-profile feeding tube made of high-quality medical silicone.

- Small, low in profile and discreet for more mobility and quality of life
- Large selection of sizes and lengths for individual and perfect fit
- Soft, flexible silicone for high patient comfort
- Better patient safety due to high-quality materials, free of harmful substances
- Apple-shaped inner retention balloon with high balloon stability for long usage, secure fixation and to avoid leakage




## The danuButton® is supplied in three different sets:

The danuButton® Premium Set and the danuButton® Compact Set include a comprehensive range of ENFit® accessories for convenient and easy replacement of the gastrostomy button. The danuButton® Compact XS is a set with reduced accessories.



| Recommended balloon filling volume |                      |
|------------------------------------|----------------------|
| Shaft size in Fr                   | Balloon volume in ml |
| 12                                 | 3.0                  |
| 14                                 | 4.0                  |
| 16                                 | 5.0                  |
| 18                                 | 5.0                  |
| 20                                 | 8.0                  |
| 24                                 | 8.0                  |



## The danumed® Extension Set

The Extension Set is the connecting piece between the danuButton® and usually an enteral syringe for the administration of enteral nutrition, liquids or medication.

The Extension Sets included in the danuButton® Premium and Compact Sets are made of kink-resistant, biocompatible polyurethane. They are equipped with ENFit® connectors according to DIN EN ISO 80369-3 and danumed's innovative rotating ring, the ENSwivel®.



## Set content of the danuButton® Premium Set:

- 1 danuButton® Gastrostomy Button
- 2 danumed® Extension Set for Gastrostomy Button, 30 cm, Bolus-ENFit®, ENSwivel®, straight button connector
- 3 danumed® Extension Set for Gastrostomy Button, 30 cm, Y-ENFit®, ENSwivel®, angled button connector
- 4 danumed® Introducer
- 5 danumed® Balloon Syringe 10 ml
- 6 danumed® Enteral Syringe, 60 ml, ENFit®
- 7 optiLube Lubricant
- 8 4 Gauze pads



## danuButton® Premium Set Product Codes:

| Shaft size 12 French |              |
|----------------------|--------------|
| Length in cm         | REF          |
| 0.8                  | DA 1041 1208 |
| 1.0                  | DA 1041 1210 |
| 1.2                  | DA 1041 1212 |
| 1.5                  | DA 1041 1215 |
| 1.7                  | DA 1041 1217 |
| 2.0                  | DA 1041 1220 |
| 2.3                  | DA 1041 1223 |
| 2.5                  | DA 1041 1225 |
| 2.7                  | DA 1041 1227 |
| 3.0                  | DA 1041 1230 |
| 3.5                  | DA 1041 1235 |
| 4.0                  | DA 1041 1240 |
| 4.5                  | DA 1041 1245 |
| 5.0                  | DA 1041 1250 |

| Shaft size 14 French |              |
|----------------------|--------------|
| Length in cm         | REF          |
| 0.8                  | DA 1041 1408 |
| 1.0                  | DA 1041 1410 |
| 1.2                  | DA 1041 1412 |
| 1.5                  | DA 1041 1415 |
| 1.7                  | DA 1041 1417 |
| 2.0                  | DA 1041 1420 |
| 2.3                  | DA 1041 1423 |
| 2.5                  | DA 1041 1425 |
| 2.7                  | DA 1041 1427 |
| 3.0                  | DA 1041 1430 |
| 3.5                  | DA 1041 1435 |
| 4.0                  | DA 1041 1440 |
| 4.5                  | DA 1041 1445 |
| 5.0                  | DA 1041 1450 |
| 5.5                  | DA 1041 1455 |
| 6.0                  | DA 1041 1460 |
| 6.5                  | DA 1041 1465 |

| Shaft size 16 French |              |
|----------------------|--------------|
| Length in cm         | REF          |
| 0.8                  | DA 1041 1608 |
| 1.0                  | DA 1041 1610 |
| 1.2                  | DA 1041 1612 |
| 1.5                  | DA 1041 1615 |
| 1.7                  | DA 1041 1617 |
| 2.0                  | DA 1041 1620 |
| 2.3                  | DA 1041 1623 |
| 2.5                  | DA 1041 1625 |
| 2.7                  | DA 1041 1627 |
| 3.0                  | DA 1041 1630 |
| 3.5                  | DA 1041 1635 |
| 4.0                  | DA 1041 1640 |
| 4.5                  | DA 1041 1645 |
| 5.0                  | DA 1041 1650 |

| Shaft size 18 French |              |
|----------------------|--------------|
| Length in cm         | REF          |
| 0.8                  | DA 1041 1808 |
| 1.0                  | DA 1041 1810 |
| 1.2                  | DA 1041 1812 |
| 1.5                  | DA 1041 1815 |
| 1.7                  | DA 1041 1817 |
| 2.0                  | DA 1041 1820 |
| 2.3                  | DA 1041 1823 |
| 2.5                  | DA 1041 1825 |
| 2.7                  | DA 1041 1827 |
| 3.0                  | DA 1041 1830 |
| 3.5                  | DA 1041 1835 |
| 4.0                  | DA 1041 1840 |
| 4.5                  | DA 1041 1845 |
| 5.0                  | DA 1041 1850 |
| 6.0                  | DA 1041 1860 |

| Shaft size 20 French |              |
|----------------------|--------------|
| Length in cm         | REF          |
| 0.8                  | DA 1041 2008 |
| 1.0                  | DA 1041 2010 |
| 1.2                  | DA 1041 2012 |
| 1.5                  | DA 1041 2015 |
| 1.7                  | DA 1041 2017 |
| 2.0                  | DA 1041 2020 |
| 2.3                  | DA 1041 2023 |
| 2.5                  | DA 1041 2025 |
| 2.7                  | DA 1041 2027 |
| 3.0                  | DA 1041 2030 |
| 3.5                  | DA 1041 2035 |
| 4.0                  | DA 1041 2040 |
| 4.5                  | DA 1041 2045 |
| 5.0                  | DA 1041 2050 |

| Shaft size 24 French |              |
|----------------------|--------------|
| Length in cm         | REF          |
| 0.8                  | DA 1041 2408 |
| 1.0                  | DA 1041 2410 |
| 1.2                  | DA 1041 2412 |
| 1.5                  | DA 1041 2415 |
| 1.7                  | DA 1041 2417 |
| 2.0                  | DA 1041 2420 |
| 2.3                  | DA 1041 2423 |
| 2.5                  | DA 1041 2425 |
| 2.7                  | DA 1041 2427 |
| 3.0                  | DA 1041 2430 |
| 3.5                  | DA 1041 2435 |
| 4.0                  | DA 1041 2440 |
| 4.5                  | DA 1041 2445 |
| 5.0                  | DA 1041 2450 |



## **danumed Medizintechnik GmbH**

*Individual solutions are our standard.*

We are a family-run, innovative, German medical technology company with high quality standards, many years of expertise in enteral nutrition and a clear customer focus.

Individual customer requests are met with flexibility and reliability.  
Personal, competent, collaborative.

danumed Medizintechnik GmbH • Stadthof 18 • 93059 Regensburg • Germany  
Tel + 49 (0) 941/ 890 588-80 • Fax + 49 (0) 941/ 890 588-81

[info@danumed.com](mailto:info@danumed.com) • [www.danumed.com](http://www.danumed.com)