



danumed® Enteral Syringes, Single Use

danumed® Enteral Syringes, Single Use, with female ENFit® connector according to DIN EN ISO 80369-3

For the administration of enteral feeds, medicines and liquids, for decompression and flushing

Extensive range of sizes for safe use in neonates, children and adults



COMFORT

Comfortable handling for the user

- **Transparent syringe cylinder** for viewing the content
- **Precisely legible** and wipe-proof scaling
- **Secure hold** with ribbed grip plate
- **Smooth**, even plunger movement
- **Noticeable safety stop** to prevent unintentional removal of the plunger from the syringe cylinder

danumed® Syringe Cap

One size fits all:
Universal syringe cap for all sizes and versions of danumed® Enteral Syringes.



danumed® Straw

Round-tip enteral straw for safe and clean withdrawal of enteral nutrition, medication and liquids.

With ENFIT® connector.

Available in different sizes and lengths.



SAFETY

Excellent safety for the patient

- **ENFit® connector** to reduce risk of tubing misconnections
- **"Low-Dose" Tip** in 1 ml, 3 ml and 5 ml syringes
- High quality biocompatible material **free of latex, PVC, DEHP and BPA**
- **Validated grading** for high measurement accuracy



danumed® Enteral Syringes, Single Use, are individually sterile packed and available in standard sizes of 1 ml, 3 ml, 5 ml, 10 ml, 20 ml, 60 ml and 100 ml. Further sizes are available on request.

Safe products according to the standards:
DIN EN ISO 80369-3, DIN EN ISO 7886-1,
DIN EN ISO 20695

danumed® Enteral Syringes, Single Use

| REF | Product description | Size | Pieces/ sales unit | Latex- free | PVC- free | DEHP- free | Bisphenol- A-free |
|--------------|---|--------|-----------------------|----------------|--------------|---------------|----------------------|
| DA 1021 0110 | danumed® Enteral Syringe, single use, with centric tip, low-dose tip and ENFit® | 1 ml | 100 | ● | ● | ● | ● |
| DA 1021 0310 | danumed® Enteral Syringe, single use, with centric tip, low-dose tip and ENFit® | 3 ml | 100 | ● | ● | ● | ● |
| DA 1021 0510 | danumed® Enteral Syringe, single use, with centric tip, low-dose tip and ENFit® | 5 ml | 100 | ● | ● | ● | ● |
| DA 1021 0500 | danumed® Enteral Syringe, single use, with centric tip and ENFit® | 5 ml | 100 | ● | ● | ● | ● |
| DA 1021 1000 | danumed® Enteral Syringe, single use, with centric tip and ENFit® | 10 ml | 100 | ● | ● | ● | ● |
| DA 1021 2000 | danumed® Enteral Syringe, single use, with centric tip and ENFit® | 20 ml | 100 | ● | ● | ● | ● |
| DA 1021 6001 | danumed® Enteral Syringe, single use, with eccentric tip and ENFit® | 60 ml | 100 | ● | ● | ● | ● |
| DA 1021 9900 | danumed® Enteral Syringe, single use, with centric tip and ENFit® | 100 ml | 25 | ● | ● | ● | ● |

All products are disposable and individually sterile packed. Further sizes are available on request.

danumed® Syringe Caps and Enteral Straws

| REF | Product description | Size/ length | Pieces/ sales unit | Latex- free | PVC- free | DEHP- free | Bisphenol- A-free |
|--------------|--|------------------|-----------------------|----------------|--------------|---------------|----------------------|
| DA 1021 0055 | danumed® Syringe Caps for Enteral Syringes, single use, ENFit® | one size | 100 | ● | ● | ● | ● |
| DA 1021 1210 | danumed® Enteral Straw for milk, ENFit® | 12 Fr/ 100 mm | 30 | ● | ● | ● | ● |
| DA 1021 1215 | danumed® Enteral Straw for milk, ENFit® | 12 Fr/ 150 mm | 30 | ● | ● | ● | ● |
| DA 1021 1604 | danumed® Enteral Straw for milk, ENFit® | 16 Fr/ 40 mm | 30 | ● | ● | ● | ● |
| DA 1021 1615 | danumed® Enteral Straw for milk, ENFit® | 16 Fr/ 150 mm | 30 | ● | ● | ● | ● |
| DA 1021 0605 | danumed® Enteral Straw for medication, ENFit® | 6 Fr/ 50 mm | 30 | ● | ● | ● | ● |

All products are disposable and individually sterile packed.

Related products to danumed® Enteral Syringes, Single Use:



danumed®
Syringe Cap



danumed®
Straw



danumed® Extension Set
for Gastrostomy Button



danuButton®
Gastrostomy Button

danumed Medizintechnik GmbH –
Individual solutions are our standard.

We are an innovative, German medical technology company with the highest quality standards, many years of expertise in enteral nutrition and clear customer focus. We also meet individual customer requests with flexibility and reliability. Personal, competent, collaborative.

danumed Medizintechnik GmbH • Stadthof 18 • 93059 Regensburg • Germany
Tel + 49 (0) 941/ 890 588-80 • Fax + 49 (0) 941/ 890 588-81

info@danumed.com • www.danumed.com



INNOVATORS IN ENTERAL CARE