



# danumed® Gravity FeedBag

Gravity Feeding Set – Enteral Bag 1500 ml with firmly connected tube system

For the administration of enteral nutrition, liquids and medication

## COMFORT

- **Large bag volume** of 1500 ml
- With large topfiller for **convenient filling** of the bag
- **Large tube lumen** for good flow rate
- **Three-way stopcock** for easy and comfortable administration of medication
- Transparent **drip chamber**
- **Roller clamp** for flow rate adjustment
- **Inscribable bag material** for clear patient identification
- **Scaling** in 100ml-steps

## SAFETY

- **High quality material** free of latex, DEHP and BPA
- **ENFit®** connectors in accordance with DIN EN ISO 80369-3 to reduce the risk of tubing misconnections and enhance patient safety
- **Protective caps** at the ENFit® connectors



REF	Product description	Pieces/sales unit
DA 1051 0210	danumed® GravityFeedBag	30 pieces

Each danumed® Gravity FeedBag is individually sterile packed. 30 pieces in a shipping box.

danumed® Gravity FeedBags are  
 • Latex-free • DEHP-free • BPA-free